

HospiCare 70

Isopropyl Alcohol 70% W Hospital Grade Disinfectant Wipes



Hard Surface Disinfection

100 sheets (20 x 18 cm)

Specially formulated with 70% Isopropyl Alcohol to effectively kill 99.99% harmful germs. Your disinfection essential for hard surfaces.

Effective against COVID-19

TGA APPROVED - AUST L: 364350









Common FAQs:

Q: What are the main features of HospiCare 70 Alcohol Wipes?

A: HospiCare 70 Alcohol Wipes are specially formulated with 70% Isopropyl Alcohol to effectively kill 99.99% harmful germs including COVID-19. Your disinfection essential for hard surfaces. It comes with a convenient and resealable lid for easy dispensing and to prevent the wipes from drying out.

Q: When and where can we use the alcohol wipes?

A: These HospiCare 70 wipes are ideal for disinfecting hard surfaces like worktops, railings, lifts, door handles, etc. It is a perfect solution for hard surface applications in hospitals, nursing homes, child-care, veterinary clinics, dental surgeries, laboraties and healthcare facilities.

Instruction for Use:



Step 1: Wipe surface in an S shape. Let the surface air dry.



Step 2: Use a new wipe for a new surface.



Step 3: Discard used wipes in a waste bin.

Q: What kind of material is used?

A: These wipes are made from a spunlace material, which is tear-resistant and soft to touch. By utilising this medium in our wipes, we increase moisture retention levels assisting in the disinfecting of hard surfaces effectively and efficiently.

Q: What is the packing size and paper size?

A: There are 100 sheets per packet and the wipe size is 20 x 18 cm.

Q: What is the shelf life of the wipes?

A: The shelf life is 2 years from date of manufacturing.

Q: How to prevent the wipes from drying?

A: Always remember to reseal the packet immediately after removing a wipe. Do not place in direct sunlight or close to any heat source. Store these wipes in a cool, dry location below 30C.





Carton External Size: 358 X 233 X 312mm